

# NGEE-Tropics ENSO A<sub>c</sub> protocol for Red Blue light source

*This protocol is for measuring A-c<sub>i</sub> curves with the LiCOR 6400XT and the 2x3cm Red Blue light source, software version 6.3.2, adapted from previous NGEA-Arctic protocols by Kim Ely and Alistair Rogers, February 2016 (updated May 2016).*

## Instruments & Equipment

### Gas Exchange

- LI-COR 6400XT gas exchange system with a R/B light source
- Sufficient batteries and chargers for the 6400 to allow measurements to be made without interruptions for charging.
- Gaskets: 2x3 white (2 day<sup>-1</sup> machine<sup>-1</sup>) and black (1 day<sup>-1</sup> machine<sup>-1</sup>), 3-hole (1 per campaign<sup>-1</sup> machine<sup>-1</sup>)
- CO<sub>2</sub> cartridges (2 day<sup>-1</sup> machine<sup>-1</sup>)
- Chemicals: soda lime, drierite
- Zero gas (UHP nitrogen), regulator, tubing
- Spare parts (see list)
- Bar codes
- Sample envelopes
- Medium ziplocks (4 x 6 inch) for holding barcoded envelop on machine
- Laptop or other alternative for downloading data
- Download ethernet adaptor card

### Advance Preparation

1. Ensure instrument is recently calibrated (<2 years)
2. Ensure most recent software is loaded (6.3.2)
3. Load NGEA Tropics config files
4. Ensure chemicals and CO<sub>2</sub> cartridges are in place at field site
5. Ensure any safety training is complete e.g. working at height
6. Instruments zero with zero gas at start of campaign

### Sampling Plan

Branches will be collected predawn / very early in the day. After removal from the tree branches will be placed in a tall bucket (kitchen trash can). A second cut will be made >1 m from the first to remove embolism. Using a transfer cup the branches will be transferred to buckets and kept in the shade for A<sub>c</sub>i measurements.

## **GAS EXCHANGE- ACi**

### **Start of Day**

1. Turn on (CO<sub>2</sub> on Scrub, Desiccant on Bypass)
2. Open config NGEE\_Tropics\_RB\_ACi
3. Follow warm up procedures
4. Open logfile using the following format
  - a. date: YYYYMMDD
  - b. site: PNM or SL
  - c. measurement type Aci
  - d. instrument ID: i.e. Derek, Mariano, Andy, Jorge, Bernie, Johnnye.g. 20160217\_SL\_ACi\_Bernie

### **Warm up procedures**

#### ***Immediately***

1. CO<sub>2</sub> cylinder in **\*replace after 6h - set a timer\***
2. Temperature (block, air, leaf) all within 2 °C
3. PARout - responds to light
4. PARin - light comes on and is stable
5. Pressure - OK (c. 100 kPa at sea level)
6. Leaf fan working
7. Flowmax > 750 μmol s<sup>-1</sup>
8. Flow restrictions - drops <15 μmol s<sup>-1</sup> on full scrub
9. Set Flowzero (pump off, leaf fan off)

#### ***After >10 minutes (full scrub desiccant and soda lime)***

10. Check CO<sub>2</sub> zero (± 5 μmol mol<sup>-1</sup> of zero)
11. Check H<sub>2</sub>O zero (0.2 mmol mol<sup>-1</sup> and falling after 2 minutes)
12. Tleaf - responds
13. Tleaf zero (block and leaf within 0.1 °C)
14. Leak test - (CO<sub>2</sub>R to 400 mol mol<sup>-1</sup>, Flow to 200 μmol s<sup>-1</sup>) blow through tube around all gaskets and seals - CO<sub>2</sub>S should change by <1 μmol mol<sup>-1</sup>
15. Match

### **Logging and sample tracking**

1. The NGEE Tropics config will prompt you to enter information in the following fields when recording a manual log (new measurements menu level 1, f1)
  - a. barcode: e.g. BNL10854 (all BNL bar codes have human readable five digit numbers, proceed this number with the BNL prefix)

- b. species: there are pull down menus for species IDs e.g. CASTEL
  - c. location: e.g. San Lorenzo
  - d. machine name: e.g. Bernie
  - e. serial number: e.g. PSC-3613
2. Keep baggie with barcode on machine during measurement
3. After each measurement write the barcode on the leaf with a sharpie or tie flagging tape marked with barcode for small leaves (to not compromise later measurements) - leaf will be collected for spectral measurement and tissue sampling later.
4. Record any notes or comments in log file using "Add Remark" function on Licor.

### Instrument stability

1. Set chamber conditions: flow  $500 \mu\text{mol s}^{-1}$ , just saturating PAR ( $2000 \mu\text{mol m}^{-2} \text{s}^{-1}$ ),  $\text{CO}_2\text{R} = 400 \text{ ppm}$ ,  $T_{\text{block}}$  at  $T_{\text{air}} + 1$  or  $2^\circ \text{C}$  (dew point buffer).
2. Wait for instrument stability
3. Match

### Leaf acclimation

1. Insert leaf and close chamber.
2. Leak test with a straw ( $<1 \mu\text{ml}$  departure from set point in  $\text{CO}_2\text{S}$ )
3. If necessary fix leaks by reclamping, tightening, removing tension on leaf or using Molycote 111 on gaskets
4. If leaf does not fill chamber adjust Area with an estimate.
5. Monitor A and  $g_s$  for a **\*minimum of 20 minutes\*** and then ensure stability (flat lines for 5 minutes) - some species took  $>45$  minutes to stabilize.
  - a. Check dew points -  $T_{\text{block}}$  should be  $1\text{-}2^\circ\text{C}$  above the dew points for sample and reference (see graph F)
  - b. Check VPD  $< 1.5 \text{ kPa}$
  - c. Check  $\Delta\text{CO}_2 > 10 \mu\text{mol s}^{-1}$  (adjust flow if necessary  $500\text{-}200 \mu\text{mol s}^{-1}$ )
  - d. Match after any adjustments

### Autoprogram

1. Record a manual log (see Logging and sample tracking) before launching the Autoprogram to load up the metadata that will be added to the data from the auto program
2. When leaf is stable, launch autoprogram "TEST\_Arctic\_Aci\_2015" (I know! - didn't have the chance to rename) Note on Johnny launch "A-Ci Curve2". (New Measurements menu 5)
3. Append to current log? Yes.
4. Start.

5. Monitor progress of A-c<sub>i</sub> curve. View graphs and “Clear all” to only see log points on the current curve.
6. After program has finished (\* next to Autoprogram menu disappears) mark up leaf if area correction is required then release from chamber. Flag leaf for collection for spectra measurement.

*Expected values for A<sub>Ci</sub> (based on Feb diurnal measurements)*

<b>Code</b>	<b>Species</b>	<b>Photosynthesis</b>	<b>Conductance</b>
CASTEL	Castilla elastica	12	0.1
LUEHSE	Luehea seemannii	22	0.3
ANACEX	Anacardium excelsum	10	0.08
CORDAL	Cordia alliodora	10	0.1
CALYCA	Calycophyllum candidissimum	10	0.1
FICUIN	Ficus insipida	18	0.2
ALIBED	Pseudosamanea guachapele	20	0.25
ANTITR	Pittoniotis trichantha	16	0.25
PSE1SE	Pseudobombax septenatum		
TERMAM	Terminalia amazonia	18	0.17
TOCOPI	Tocoyena pittieri	9	0.1
CARAGU	Carapa guianensis	12	0.15
TACHVE	Tachigali versicolor	15	0.5
VOCHFE	Vochysia ferruginea	20	0.45
VIROSP	Virola multiflora	7	0.05
MICOBO	Miconia borealis	15	0.25
APEIME	Apeiba membranacea	20	0.5
GUATDU	Guatteria dumetorum	15	0.6